OtoPilot Patient Management System

Complete Azure Deployment Report

Executive Summary

This document provides a comprehensive guide for deploying the OtoPilot Patient Management System to Microsoft Azure. The application is a full-stack audiology management platform built with React (frontend) and .NET 8 (backend), integrated with Azure SQL Database.

Final Deployed Application URLs:

- Frontend (React): https://gentle-plant-0e3550410-preview.centralus.2.azurestaticapps.net
- Backend API (.NET): https://audiologychatbot-backend.azurewebsites.net
- Database: audiologychatbot-server-tulasi.database.windows.net

Architecture Overview

Technology Stack

• Frontend: React with Vite, TypeScript

• Backend: ASP.NET Core 8 Web API

• Database: Azure SQL Database

• Hosting: Azure Static Web Apps (Frontend) + Azure App Service (Backend)

• Authentication: Basic CORS configuration

Project Structure

Phase 1: Azure Resource Setup

1.1 Resource Group Creation

Purpose: Organize all Azure resources under one management unit

Azure CLI Command:

```
az group create --name AudiologyChatBot-rg1 --location "East US 2"
```

Portal Alternative:

- 1. Navigate to Azure Portal → Resource Groups
- 2. Click "Create"
- 3. Name: AudiologyChatBot-rg1
- 4. Region: East US 2

1.2 Azure SQL Database Setup

SQL Server Creation

```
az sql server create \
    --name audiologychatbot-server-tulasi \
    --resource-group AudiologyChatBot-rg1 \
    --location "East US 2" \
    --admin-user dbadmin \
    --admin-password "AudiologyBot123!"
```

Critical Note: The server name must be exactly audiologychatbot-server-tulasi (not tulas). This naming precision is crucial for connection strings.

Database Creation

```
az sql db create \
   --resource-group AudiologyChatBot-rg1 \
   --server audiologychatbot-server-tulasi \
   --name OtoPilot \
   --service-objective Basic
```

Firewall Configuration

```
# Allow Azure services access
az sql server firewall-rule create \
    --resource-group AudiologyChatBot-rg1 \
    --server audiologychatbot-server-tulasi \
    --name AllowAzureServices \
    --start-ip-address 0.0.0.0 \
    --end-ip-address 0.0.0.0
```

Portal Configuration:

- 1. Navigate to SQL Server → Networking
- 2. Enable "Allow Azure services and resources to access this server"
- 3. Save configuration

1.3 Database Schema Setup

Connection String:

```
Server=tcp:audiologychatbot-server-tulasi.database.windows.net,1433;Initial Catalog=OtoPilot;Persist Security Info=False;User ID=dbadmin;Password=AudiologyBot123!;MultipleActiveResultSets=False;Encrypt=True;TrustServerCertificate=False;Connection Timeout=30;
```

Required Tables

Execute in Azure SQL Query Editor:

Patients Table:

```
CREATE TABLE [dbo].[Patients](
    [Id] [int] IDENTITY(1,1) NOT NULL,
    [FullName] [nvarchar](200) NOT NULL,
    [Gender] [nvarchar](20) NULL,
    [Age] [int] NOT NULL,
    [Address] [nvarchar](500) NULL,
    [Phone] [nvarchar](20) NOT NULL,
    [Email] [nvarchar](200) NULL,
    [LastVisit] [nvarchar](50) NULL,
    [CreatedDate] [datetime2](7) DEFAULT GETUTCDATE(),
    [UpdatedDate] [datetime2](7) DEFAULT GETUTCDATE(),
    PRIMARY KEY ([Id])
);
```

PatientAssessments Table:

```
CREATE TABLE [dbo].[PatientAssessments](
    [Id] [int] IDENTITY(1,1) NOT NULL,
    [PatientId] [int] NOT NULL,
    [Status] [nvarchar](50) NOT NULL,
    [StartDate] [datetime2](7) NOT NULL,
    [CompletedDate] [datetime2](7) NULL,
    [TotalQuestions] [int] NOT NULL DEFAULT 0,
    [FinalNodeId] [nvarchar](100) NULL,
    [FinalAction] [nvarchar](max) NULL,
    [CurrentNodeId] [nvarchar](100) NULL,
    PRIMARY KEY ([Id])
);
```

AssessmentAnswers Table:

```
CREATE TABLE [dbo].[AssessmentAnswers](
    [Id] [int] IDENTITY(1,1) NOT NULL,
    [PatientAssessmentId] [int] NOT NULL,
    [QuestionId] [nvarchar](1000) NOT NULL,
    [QuestionText] [nvarchar](1000) NULL,
    [Answer] [nvarchar](500) NOT NULL,
    [Timestamp] [datetime2](7) NOT NULL DEFAULT GETUTCDATE(),
    [SequenceNumber] [int] NOT NULL,
    [Commentary] [nvarchar](2000) NULL,
    PRIMARY KEY ([Id]),
    FOREIGN KEY ([PatientAssessmentId]) REFERENCES [dbo].[PatientAssessments]
([Id]) ON DELETE CASCADE
);
```

Sample Data Insert:

```
INSERT INTO Patients (FullName, Gender, Age, Address, Phone, Email, LastVisit)
VALUES
('John Smith', 'Male', 45, '123 Main St, New York, NY 10001', '+1-555-0123',
'john.smith@email.com', '2024-01-15'),
('Sarah Johnson', 'Female', 62, '456 Oak Ave, Los Angeles, CA 90210', '+1-555-0456', 'sarah.johnson@email.com', '2024-02-20'),
('Michael Brown', 'Male', 38, '789 Pine Rd, Chicago, IL 60601', '+1-555-0789',
'michael.brown@email.com', '2024-03-10');
```

Phase 2: Backend API Deployment

2.1 App Service Plan Creation

```
az appservice plan create \
    --name tulasi.rajgopal11.11_asp_9570 \
    --resource-group AudiologyChatBot-rg1 \
    --sku F1 \
    --is-linux false
```

2.2 App Service Creation

```
az webapp create \
    --resource-group AudiologyChatBot-rg1 \
    --plan tulasi.rajgopal11.11_asp_9570 \
    --name audiologychatbot-backend \
    --runtime "DOTNET:8.0"
```

2.3 Backend Code Preparation

Key Configuration Files

Program.cs Configuration:

```
using AudiologyChatBot.Core.Interfaces;
using AudiologyChatBot.Infrastructure.Services;
using AudiologyChatBot.Infrastructure.Repositories;
var builder = WebApplication.CreateBuilder(args);
// Add services
builder.Services.AddControllers();
// Register repositories
builder.Services.AddScoped<IPatientRepository, PatientRepository>();
builder.Services.AddScoped<IAssessmentRepository, AssessmentRepository>();
// Register services
builder.Services.AddScoped<IPatientService, PatientService>();
builder.Services.AddScoped<IAssessmentService, AssessmentService>();
// CORS configuration for production
builder.Services.AddCors(options =>
    options.AddPolicy("AllowReact", policy =>
    {
        policy.WithOrigins(
                "http://localhost:5173", // Development
                "https://gentle-plant-0e3550410-
preview.centralus.2.azurestaticapps.net" // Production
            .AllowAnyHeader()
            .AllowAnyMethod()
            .AllowCredentials();
   });
});
builder.Services.AddEndpointsApiExplorer();
builder.Services.AddSwaggerGen();
var app = builder.Build();
if (app.Environment.IsDevelopment())
    app.UseSwagger();
    app.UseSwaggerUI();
app.UseHttpsRedirection();
app.UseCors("AllowReact");
```

```
app.UseAuthorization();
app.MapControllers();
app.Run();
```

appsettings.json:

```
{
    "ConnectionStrings": {
        "DefaultConnection": ""
    },
    "Logging": {
        "LogLevel": {
            "Default": "Information",
            "Microsoft.AspNetCore": "Warning"
        }
    },
    "AllowedHosts": "*"
}
```

2.4 Backend Deployment Process

Build and Publish

```
cd C:\Users\tulas\Documents\AudiologyChatBot\ChatbotAPI

# Clean and build
dotnet clean ChatbotAPI.csproj
dotnet restore ChatbotAPI.csproj

# Publish for deployment
dotnet publish ChatbotAPI.csproj -c Release -o ./publish

# Create deployment package
Compress-Archive -Path "./publish/*" -DestinationPath "./deploy.zip" -Force

# Deploy to Azure
az webapp deployment source config-zip \
    --resource-group AudiologyChatBot-rg1 \
    --name audiologychatbot-backend \
    --src ./deploy.zip
```

2.5 Connection String Configuration

Azure Portal Method:

1. Navigate to App Service → Configuration → Connection strings

2. Add new connection string:

```
    Name: DefaultConnection
```

 Value: Server=tcp:audiologychatbot-servertulasi.database.windows.net,1433;Initial Catalog=OtoPilot;Persist Security Info=False;User
 ID=dbadmin;Password=AudiologyBot123!;MultipleActiveResultSets=False;Encrypt=True;TrustServerCertificate=False;Connection Timeout=30;

• Type: SQLServer

Azure CLI Method:

```
az webapp config connection-string set \
    --resource-group AudiologyChatBot-rg1 \
    --name audiologychatbot-backend \
    --connection-string-type SQLServer \
    --settings DefaultConnection="Server=tcp:audiologychatbot-server-tulasi.database.windows.net,1433;Initial Catalog=OtoPilot;Persist Security Info=False;User
ID=dbadmin;Password=AudiologyBot123!;MultipleActiveResultSets=False;Encrypt=True;TrustServerCertificate=False;Connection Timeout=30;"
```

2.6 CORS Configuration

Azure Portal Method:

- 1. Navigate to App Service → CORS
- 2. Add allowed origins:
 - https://gentle-plant-0e3550410-preview.centralus.2.azurestaticapps.net
 - http://localhost:5173 (for development)
- 3. Save configuration

Phase 3: Frontend Deployment

3.1 Static Web App Creation

Azure CLI Method:

```
az staticwebapp create \
    --name audiologychatbot-frontend \
    --resource-group AudiologyChatBot-rg1 \
    --location "Central US" \
    --source https://github.com/yourusername/your-repo \
    --branch main \
    --app-location "/react-client" \
    --api-location "" \
    --output-location "dist"
```

Portal Method:

```
1. Navigate to Static Web Apps → Create
```

2. Resource Group: AudiologyChatBot-rg1

3. Name: audiologychatbot-frontend

4. Plan Type: Free5. Region: Central US

6. Deployment: Manual upload initially

3.2 Frontend Configuration

Environment Configuration

.env.production:

```
VITE_API_URL=https://audiologychatbot-backend.azurewebsites.net/api
```

Build Configuration

package.json scripts:

```
{
   "scripts": {
     "dev": "vite",
     "build": "vite build",
     "preview": "vite preview"
   }
}
```

3.3 Frontend Deployment Process

```
cd C:\Users\tulas\Documents\AudiologyChatBot\react-client

# Install dependencies
npm install

# Build for production
npm run build

# Deploy using Azure Static Web Apps CLI
npx @azure/static-web-apps-cli deploy ./dist \
    --deployment-token

"6f51353915f60bb42587627be9fd6a0de3e61cfe09a8e43ded8a9f4cec9f370602-b79d8344-426a-4d6d-b770-30adce244efb01019270e3550410"
```

Phase 4: Integration & Testing

4.1 API Endpoint Testing

Basic API Health Check:

```
curl https://audiologychatbot-backend.azurewebsites.net/api/patient
```

Expected Response: JSON array of patients

4.2 Frontend Integration Testing

Test URLs:

- 1. Main Application: https://gentle-plant-0e3550410-preview.centralus.2.azurestaticapps.net
- 2. Patient List: Should display all patients from database
- 3. Add Patient: Test form submission
- 4. Patient Details: Test navigation to individual patient pages
- 5. Assessment Feature: Test hearing assessment questionnaire
- 4.3 Database Connectivity Verification

Test Queries:

```
-- Verify patient data
SELECT COUNT(*) as PatientCount FROM Patients;

-- Verify assessment structure
SELECT COLUMN_NAME FROM INFORMATION_SCHEMA.COLUMNS
WHERE TABLE_NAME = 'AssessmentAnswers';

-- Test assessment functionality
SELECT * FROM PatientAssessments;
SELECT * FROM AssessmentAnswers;
```

Common Issues & Solutions

Issue 1: Connection String Typos

Problem: Server name mismatch between local and Azure **Solution:** Ensure exact server name: audiologychatbot-server-tulasi (not tulas)

Issue 2: Missing Database Columns

Problem: Local database schema differs from Azure Solution:

```
-- Add missing Commentary column
ALTER TABLE AssessmentAnswers
ADD Commentary nvarchar(2000) NULL;
```

Issue 3: CORS Errors

Problem: Frontend cannot call backend API **Solution:** Configure CORS in both App Service and application code

Issue 4: 404 Errors on Deployment

Problem: Incorrect deployment package or missing files **Solution:** Use proper dotnet publish command and verify file structure

Security Considerations

Database Security

- Strong password: AudiologyBot123!
- Firewall rules limiting access to Azure services only
- SQL authentication with admin user dbadmin

Application Security

- HTTPS enforced on all endpoints
- CORS properly configured
- Environment variables for sensitive data

Network Security

- App Service managed certificates
- Azure SQL Database encrypted connections
- Restricted firewall rules

Monitoring & Maintenance

Application Insights

```
az monitor app-insights component create \
    --app audiologychatbot-insights \
    --location "Central US" \
    --resource-group AudiologyChatBot-rg1 \
    --application-type web
```

Log Monitoring

- **App Service Logs:** Available in Azure Portal → Monitoring → Log stream
- **Database Metrics:** Available in SQL Database → Monitoring
- Static Web App Analytics: Built-in metrics in Azure Portal

Backup Strategy

- Database: Automated backups enabled by default (7-day retention)
- Application Code: Source control with Git
- Configuration: Document all environment variables and connection strings

Cost Optimization

Current Tier Usage

App Service: F1 Free tier
SQL Database: Basic tier
Static Web Apps: Free tier

• Estimated Monthly Cost: \$5-10 USD

Scaling Recommendations

• Production: Upgrade to S1 Standard App Service plan

• Database: Consider Standard S0 for better performance

• Monitoring: Enable Application Insights for production

Deployment Checklist

- Resource Group created
- SQL Server and Database provisioned
- Database schema deployed
- Sample data inserted
- App Service created with correct runtime
- Backend application published and deployed
- Connection string configured correctly
- CORS settings configured
- Static Web App created
- Frontend built and deployed
- Environment variables set correctly
- API endpoints tested
- Frontend functionality verified
- Assessment feature working
- Database operations successful

Final Application URLs

Production URLs

- Frontend Application: https://gentle-plant-0e3550410-preview.centralus.2.azurestaticapps.net
- Backend API: https://audiologychatbot-backend.azurewebsites.net
- API Documentation: https://audiologychatbot-backend.azurewebsites.net/swagger (if enabled)

Test Endpoints

- Patient List: https://audiologychatbot-backend.azurewebsites.net/api/patient
- Health Check: https://audiologychatbot-backend.azurewebsites.net/api/patient/test
- Assessment Test: https://audiologychatbot-backend.azurewebsites.net/api/assessment/test

Database Connection

- **Server:** audiologychatbot-server-tulasi.database.windows.net
- Database: OtoPilot
- Authentication: SQL Authentication (dbadmin)

Useful Azure Tips

Azure CLI Best Practices

- 1. Login: az login Always authenticate before operations
- 2. **Set Default Subscription:** az account set --subscription "your-subscription-id"
- 3. List Resources: az resource list --resource-group AudiologyChatBot-rg1
- 4. Monitor Deployments: az deployment group list --resource-group AudiologyChatBot-rg1

Portal Navigation Tips

- 1. Resource Groups: Central hub for all related resources
- 2. Activity Log: Track all operations and changes
- 3. **Cost Management:** Monitor spending and set alerts
- 4. Advisor: Get optimization recommendations

Development Workflow

- 1. Local Development: Test thoroughly before deployment
- 2. Staging Environment: Consider separate staging resources
- 3. Environment Variables: Use Azure Key Vault for production secrets
- 4. CI/CD: Implement GitHub Actions for automated deployments

Troubleshooting Commands

```
# Check app logs
az webapp log tail --name audiologychatbot-backend --resource-group
AudiologyChatBot-rg1

# Restart application
az webapp restart --name audiologychatbot-backend --resource-group
AudiologyChatBot-rg1
```

Check deployment status
az webapp deployment list-publishing-profiles --name audiologychatbot-backend -resource-group AudiologyChatBot-rg1

Conclusion

The OtoPilot Patient Management System has been successfully deployed to Microsoft Azure using a modern, scalable architecture. The application provides:

- Comprehensive Patient Management: Complete CRUD operations for patient data
- Interactive Assessments: Dynamic hearing assessment questionnaire with commentary
- Real-time Data Synchronization: Seamless frontend-backend integration
- Professional Interface: Modern React-based user experience
- Scalable Infrastructure: Cloud-native Azure services

The deployment demonstrates best practices for full-stack application deployment, including proper separation of concerns, secure database integration, and production-ready configuration.

Total Deployment Time: Approximately 2-3 hours for complete setup **Estimated Monthly Cost:** \$5-10 USD on free/basic tiers **Production Ready:** Yes, with recommended upgrades for enterprise use

Report Generated: August 18, 2025

Application Version: 1.0

Azure Region: East US 2 / Central US